

**TRIOLOGY[®]
MULTIFOCALS**



Superior Optics



TRIOLOGY[®]
LENSES

**Impact Resistant
&
Chemical Resistant**

**Trilogy[®]
Flat-Top 28
Bifocals**



**NEW:
Trilogy[®]
Flat-Top 35
Bifocals**

Lightweight

**NEW:
Trilogy[®]
7 x 28
Trifocals**

**YOUNGER
OPTICS** 
The Optical Lens Innovators

LAB ANNOUNCEMENT:

Trilogy® Multi-Focals now available in FT35 & 7x28 from Younger Optics

Trilogy® FT28 Bifocal - Clear, Hardcoated

Ø	NOMINAL CURVE	TRUE CURVE	BACK CURVE	EDGE	CENTER	SEG DROP	SEG INSET	SAG AT 50	RADIUS	RX RANGE
76	2.25	2.26	4.25	14.5	12.5	5.0	6.0	1.33	234.280	-8.00D to +6.50D
76	4.25	4.28	6.25	14.5	11.5	5.0	6.0	2.55	123.830	
76	6.25	6.30	6.25	8.0	8.0	5.0	6.0	3.80	84.195	
76	8.25	8.34	6.25	8.5	11.5	5.0	6.0	5.12	63.573	

Add Range 1.00 D to 3.50 D

NEW Trilogy® FT35 Bifocal - Clear, Hardcoated

Ø	NOMINAL CURVE	TRUE CURVE	BACK CURVE	EDGE	CENTER	SEG DROP	SEG INSET	SAG AT 50	RADIUS	RX RANGE
76	2.25	2.26	4.25	14.5	12.5	5.0	6.0	1.33	234.280	-8.00D to +6.50D
76	4.25	4.28	6.25	14.5	11.5	5.0	6.0	2.55	123.830	
76	6.25	6.29	6.25	8.0	8.0	5.0	6.0	3.80	84.195	
76	8.25	8.34	6.25	8.5	11.5	5.0	6.0	5.12	63.573	

Add Range 1.00 D to 3.50 D

NEW Trilogy® 7 x 28 Trifocal - Clear, Hardcoated

Ø	NOMINAL CURVE	TRUE CURVE	BACK CURVE	EDGE	CENTER	SEG DROP	SEG INSET	SAG AT 50	RADIUS	RX RANGE
76	4.25	4.30	6.28	14.2	11.8	4.0	4.0	2.56	123.272	-7.00D to +6.50D
76	6.25	6.30	6.28	7.7	8.0	4.0	4.0	3.80	84.103	
76	8.25	8.36	6.28	7.8	11.3	4.0	4.0	5.14	63.397	

Add Range 1.50 D to 3.50 D

To order, call 1-800-366-5367

YOUNGER OPTICS 
The Optical Lens Innovators