

NOW AVAILABLE FROM
YOUNGER OPTICS
 TO WHOLESALE LABS

ULTRA
 HIGH INDEX

1.74
 CLEAR

MADE FROM MR-174™ MATERIAL



PREMIUM CHOICE
 for high prescriptions

THINNER LENSES
 for cosmetic appeal

RECOMMENDED
 for rimless frames

WIDE RANGE
 of base curves

AVAILABILITY

COATED AND UNCOATED SFSV | **RX RANGE: -18.00 D TO +12.00 D**

Diameter (mm)	Nominal Curve	True Curve	Back Curve	Edge Thickness	Center Thickness	Sagitta (50 mm)	Radius (mm)
65	0.5	0.27	9.4	18.5	8.5	0.16	1962.96
65	1.0	1.04	9.4	16.8	7.5	0.61	509.62
70	1.5	1.57	8.2	15.5	7.0	0.93	337.58
75	2.5	2.47	6.9	12.5	6.1	1.46	214.57
75	3.5	3.54	4.8	8.0	6.3	2.10	149.80
70	5.0	4.78	5.6	7.0	6.0	2.86	110.83
70	6.0	6.12	6.2	7.1	7.0	3.68	86.60
70	7.0	7.31	6.9	6.9	7.5	4.45	72.47
65	8.0	8.32	6.7	6.6	8.6	5.11	63.70
65	9.0	8.95	6.7	6.6	9.5	5.53	59.24
65	10.0	9.7	6.2	6.5	10.7	6.06	54.63

WHOLESALE LABS: CONTACT YOUR YOUNGER REPRESENTATIVE
 OR CALL YOUNGER OPTICS CUSTOMER SERVICE **800-366-5367**

YOUNGER OPTICS
 The Optical Lens Innovators